



Plastic Weaner

Weiyou® Plastic weaner is made of high quality ABS Non-toxic plastic. Weiyou® Plastic weaner is used to prevent calfs to lick bellybutton of each other and such the breast. As the manufacturer in China, Weiyou® has mature production lines and supply chains, and always supply best quality product, best favourable price and best service.



Product Description

Weiyou® Plastic weaner is made of high quality ABS Non-toxic plastic. Weiyou® Plastic weaner is used to prevent calfs to lick bellybutton of each other and such the breast. Weiyou® Plastic weaner is easy to adjust the screw, can use again, and avoid the mavericks such behavior, avoid to spread the disease. Various colors are available, and OEM customization is acceptable. Vaccination to prevent disease for all kinds of animals. Free samples, and small quantity can be supplied.

Specification

Item	Plastic weaner
Usage	Prevent calfs to lick bellybutton of each other
Material	High quality ABS Non-toxic plastic
MOQ	500pcs
Use For	Cattle, Horse Sheep, Goat, etc
Color	Orange, Red, Yellow,Blue and other colors
Packing	Carton

Features

Plastic weaner for cow calf

Easy to adjust the screw,can use again

Avoid the mavericks such behavior, avoid to spread the disease

Used to prevent the calfs to such the breast

Used to prevent calfs to lick bellybutton of each other





Rational construction

It is not easy to damage, the force is more uniform and stable, will be weaned to the end

Product scope

Use for calves and lambs in weaning



Product features

The main thickened material, durable, not easy to deformation









